



reddot design award
winner 2010

BIONIK^{ZEPP} Prophylaxis

BIONIK^{ZEPP}

The new **ZEPP BIONIK**
Handle corresponds with the
demand of an anatomically ad-
justed handle for prophylaxis.

The perfect shape regarding
power transmission and sensi-
tivity enables a tactile scaling
and cleaning.

The adaptation of the prac-
tice-oriented requirements
regarding communication,
hygienics and flexibility make
the **HELMUT ZEPP BIONIK**
Handle the perfect instrument
holder not only for dental diag-
nosis and prophylaxis, but also
for surgery, implantology and
microsurgery.



HELMUT ZEPP
MEDIZINTECHNIK GMBH

NEWS

03/14

PERIODONTAL PROBES MOUTH MIRRORS

INFORMATION:

For further information on the
ZEPP BIONIK Universal
Handles please send us your
inquiry or contact us online at:
www.zepp-dental.com

All handle colors of
HELMUT ZEPP BIONIK
Universal Handles:



26.193.10	26.193.08	26.193.06	26.193.04	26.193.02
26.193.09	26.193.07	26.193.05	26.193.03	26.193.01

2

ZEPP

Periodontal Probes


partially with Color Coding

The perio-diagnostics should be done regularly to have a permanent record to make sure that gingival and especially periodontal inflammations are diagnosed early.

The benefit for the patient is that either the periodontal health can be confirmed or infection is detected early. That way, suitable dental actions can take place to cure the infection or to stop the progression.

With the new color-coded periodontal probes an even easier visualisation is guaranteed during the examination.

The green band (3 mm) of the probe stays visible at the deepest point of the sulcus of all sextants teeth. Calculus or defective restoration margins cannot be detected, the tissue of the gingiva is healthy, no bleeding after (careful) probing = no treatment necessary, further preventive care. Deeper measurements (red) are indicating a gingivitis. According to the depth, further evaluation and documentation must take place according to the rules of the PSI.



Color Periodontal Probe with Color Coding

24.454.01C

CPG12


3 - 6 - 9 - 12 mm



24.454.01

CPG 12

3 - 6 - 9 - 12 mm



24.454.02

CPNG 22

2 - 4 - 6 - 8 - 10 - 12 mm



ZEPP Periodontal Probe M2.5

26.194.08

HELMUT ZEPP BIONIK Universal
Handle, PEEK, for exchangeable inserts
M2.5 x 0.5, single-ended, yellow green

Their interchangeable tips remain firmly held in place at all times, thanks to the conical necks on the latter's threaded ends.



WHO

24.454.03

CPG 11.5 (WHO)

3.5 - 5.5 - 8.5 - 11.5 mm



ZEPP Lipcare Mouth Mirror

The Lipcare Mouth Mirror / Retractor Combination is distinguished by the following features:

- optimal ergonomic working for the user
- patient-friendly, atraumatic, extensive retraction of the cheek
- Titanium Mirror, Ø 24 mm, for best visibility
- possibility of preparation with all usual cleaning methods
- sterilizable up to 273° F / 134° C

Due to the standardized thread dimensions, all mirrors available from **ZEPP** are compatible with all handles delivered by **HELMUT ZEPP** and can of course be disinfected and sterilized by all usual methods.



24.062.24 VE 5

24.062.24. VE 1

Lipcare Mouth Mirror, Titanium,
Ø 24 mm, sterilizable

26.194.03

ZEPP BIONIK Universal Handle,
PEEK, for exchangeable inserts
M2.5, single-ended, red-purple




Front Surface Mirrors

These Mirrors have their reflective coating on their front surface, yielding undistorted imaging with no disturbing double images.


Tel.: +49 (0) 74 64/98 88-0

ZEPP BIONIK Gracey Special Cures


ZEPP Gracey Cures are special Cures that have just a single working surface or cutting edge, and are thus suitable for removing concretations or dental plaque only. The inserts are exchangeable.




24.751.101G GRA 1/2 Special Curette for front teeth, slightly angled, yellow **24.751.102G**




24.751.103G GRA 3/4 Special Curette for front teeth/premolars, slightly angled, signal orange **24.751.104G**




24.751.105G GRA 5/6 Special Curette for front teeth/premolars, more sharply angled, red-purple **24.751.106G**




24.751.107G GRA 7/8 Special Curette for premolars/molars, more sharply angled, light red-magenta **24.751.108G**




24.751.109G GRA 9/10 Special Curette for premolars/molars, very sharply angled, turquoise-bright blue **24.751.110G**




24.751.111G GRA 11/12 Special Curette for use on all mesial surfaces of premolars/molars. Skewed to allow optimal placement, signal purple **24.751.112G**



24.751.113G GRA 13/14 Special Curette for use on all distal surfaces of premolars/molars. Skewed to allow optimal placement, cobalt blue **24.751.114G**



24.751.115G GRA 15/16 Special Curette for use on all mesial surfaces of premolars/molars. Skewed like 13/14, opposite working surface, signal purple **24.751.116G**



24.751.117G GRA 17/18 Special Curette for use on premolars/molars. Triple-skewed for optimal access to distal surfaces or deep pockets. Provides good access, even when mouth openings are restricted, cobalt blue **24.751.118G**


The BIONIK Handle

The adaptation of the practice-oriented requirements regarding communication, hygienics and flexibility make the **ZEPP BIONIK** Handle the perfect instrument holder not only for dental diagnosis and prophylaxis, but also for surgery, implantology and microsurgery.



HELMUT ZEPP
MEDIZINTECHNIK GMBH

www.zepf-dental.com



ZEPP BIONIK M5 Currettes | Deep Scaling

HELMUT ZEPP M5 Currettes have a 1st shaft which is about 3 mm longer. The sharpened instrument tip is shortened as to allow a special subgingival curettage for tight and deep pockets as well as narrow root surfaces. The inserts are exchangeable.



ZEPP Safety Key

The Safety Key holds the tips firmly in place. The side wings (only **24.755.02**) protect from potential injuries due to the extremely sharp cutting edges of the tips.



ZEPP BIONIK All-purpose Curettes

Unlike Gracey Curettes, **HELMUT ZEPP** All-purpose Curettes have two working surfaces or cutting edges, and thus allow using the same instrument on teeth's mesial and distal surfaces. The inserts are exchangeable.

5

Exchanging Curette and Scaler Inserts:




Start by screwing the tip (sharp end) carefully by hand into the handle. Make sure that only the handle is turned. Not the tip.



Pull in the tip with the **ZEPP** Safety Key. With the other end, proceed as before. To set the tips, it is better to use two Safety Keys.




The Safety Key holds the tips firmly in place. The side wings protect from potential injuries due to the extremely sharp cutting edges of the tips.



24.751.102LC

Columbia University, fig. CU 2R/2L, Curette for use on all front teeth, cobalt blue


24.751.102RC



24.751.113C

Columbia University, fig. CU 13/14, Curette suitable for all types of use, cobalt blue


24.751.114C




24.751.104LC

Columbia University, fig. CU 4R/4L, Curette for use on all side teeth, cobalt blue

24.751.104RC




24.751.113MC McCall, fig. MC 13S/14S, Scaler for use on interdental spaces **24.751.114MC** on front teeth and premolars, also usable for general work. Fine, sharply pointed, cobalt blue



24.751.117MC

McCall, fig. MC 17S/18S, all-purpose Curette for use on molars, cobalt blue


24.751.118MC



24.751.101L

Langer, fig. L 1/2, Curette for use on lower molars and premolars, black


24.751.102L



24.751.103L

Langer, Fig. L 3/4, Curette for use on upper molars and premolars, yellow-green


24.751.104L



24.751.105L

Langer, Fig. L 5/6, Curette for use on upper and lower front teeth, yellow

24.751.106L




24.751.123L

M 23 Scaler for use on upper and lower side teeth and supragingival dental plaque. Slim version, cobalt blue

24.751.123R

TITANIUM CURETTES ALL-PURPOSE + M5

6



Standard

M5

HELMUT ZEPF
MEDIZINTECHNIK GMBH

Tel.: +49 (0) 74 64/98 88-0

TI ZEPF BIONIK Titanium Currettes

All-purpose Currettes with exchangeable, very delicate titanium inserts from **HELMUT ZEPF**, to remove the accumulated plaque film on the implant necks. The titanium inserts are exchangeable.




TI ZEPF BIONIK M5 Titanium Currettes

HELMUT ZEPF M5 Currettes have a 1st shaft which is about 3 mm longer. To remove the accumulated plaque film on the implant necks. The titanium inserts are exchangeable.



ZEPP BIONIK Scaler


HELMUT ZEPP Sickle Scalers have ultrafine tips for supragingival use. Unlike Jacquette scalers, they have curved facial surfaces. With blunt tips, for general use. The inserts are exchangeable.



24.751.130

Scaler for use on front teeth,
yellow-green


24.751.131



24.751.107H


Hygienist, fig. H 6/7, Sickle Scaler,
with opposed tips, for front teeth and premolars, yellow-green

24.751.106H



24.751.102CI Mini-Kaplan, fig. CI 2/3, Sickle Scaler with large facial surfaces
matched to tooth contours for removing massive accumulations of dental plaque
or plastic, cement, or adhesive residues from all interdental spaces, yellow-green


24.751.103CI



24.751.135UL

Sickle Scaler, fig. U 135, with small facial surfaces matched
to tooth contours for use in all interdental spaces, yellow-green


24.751.135UR



24.751.102TA

Taylor, fig. T 2/3, medium-sized Scaler with facial surfaces
matched to tooth contours for use in all interdental spaces, yellow-green


24.751.103TA



24.751.204L

Sickle Scaler, fig. 204 S, for removing dental plaque
from interdental spaces in the molar area, yellow-green


24.751.204R



24.751.204LD

Sickle Scaler, fig. 204 SD, for removing dental plaque
from narrow interdental spaces in the premolar and frontal area, yellow-green


24.751.204RD



24.751.123AL

Sickle Curette (Molar Scaler), Fig. M 23A,
with blunt tips for removing dental plaque from interdental spaces
in the molar area. Also suitable for general use, yellow-green

24.751.123AR



24.751.104XL

Sickle Curette (Scaler), Fig. GXC 4,
with blunt tips for general use, grey

24.751.104XR

BIONIK ID-Plugs

When placing orders for individual **BIONIK** Handles or re-ordering them, please indicate the desired ID-plug **separately** according to our listing below.

Material: PEEK

-  **26.195.01** yellow
-  **26.195.02** signal orange
-  **26.195.03** red-purple
-  **26.195.04** signal purple
-  **26.195.05** lightred-magenta
-  **26.195.06** turquoise-brightblue
-  **26.195.07** cobalt blue
-  **26.195.08** yellow-green
-  **26.195.09** grey
-  **26.195.10** black



HELMUT ZEPP
MEDIZINTECHNIK GMBH
www.zepf-dental.com



red dot design award
winner 2010

Distributed by:

BIONIK ZEPF

Communication:

- 11 different, fresh base colours
- additional ID-plugs in corresponding colors

Hygienics:

- PEEK High-Tech Plastic Material
- the shape design allows best cleaning results and ideal grip

Flexibility:

- exchangeable working tips
- the **BIONIK**-Line is expanded to the fields of surgery, implantology and microsurgery

Ergonomics:

- perfect handiness
- ideal power transmission with formerly unknown sensitivity
- relaxed working as the handle can be held in different (diameter) positions

24.990.50

ZEPF BIONIK Prophylaxis-Set 'Gracey'
consisting of Gracey 5/6, 7/8, 11/12, 13/14
Scaler 204S and 1/3 Washtray

24.990.55

ZEPF BIONIK Prophylaxis-Set 'M5 Deep
Scaling' consisting of Gracey 1/2 M5, 7/8 M5,
11/12 M5, 13/14 M5, Scaler 204S and 1/3 Washtray

24.990.60

ZEPF BIONIK Prophylaxis-Set 'Universal'
consisting of Langer 1/2, 3/4, 5/6, M23
Scaler 204S and 1/3 Washtray

BIONIK



HELMUT ZEPF
MEDIZINTECHNIK GMBH